



# INFECTION PREVENTION & CONTROL



AS-5369:2023: Australian Standard™  
Reprocessing of reusable medical  
devices (RMD) and other devices in  
health and non-health related facilities

## PROTECT YOUR PATIENTS. PROTECT YOUR PRACTICE.

RMD reprocessing is a multistep process that may include cleaning, disinfection, inspection and assembly, testing, packaging and sterilisation, transport and storage with appropriate handling to render them safe for reuse

### Minimise the risk

- **Practice hand hygiene:** In accordance with the 5 Moments of Hand Hygiene. Don suitable gloves for invasive procedures.
- **Use aseptic technique:** Maintain an aseptic field to prevent contamination during semi-critical and critical procedures.
- **Routine education recommended:** Complete infection prevention and control modules on practice as required, and in line with facility or other requirements.
- **Risk assessment & documentation:** Access and document protocols for handling and storing RMDs

### Reprocessing the transducer

#### Step 1: Cleaning the transducer

- **Dedicated reprocessing area:** Set up a clearly defined space with unidirectional workflow to prevent cross contamination
- **Approved cleaning products:** Use Class I cleaning agents listed on the ARTG, or cleaning agents listed on the NZ Medsafe website, for safe and effective cleaning
- **Remove all debris from the transducer:** Including gel, organic matter, and residual cleaning agent to achieve a clean and dry transducer prior to HLD

#### Step 2: High Level Disinfection (HLD)

- **Validated disinfection:** Only use Class IIb disinfection products listed on the ARTG, or disinfection products listed on the NZ Medsafe website, to meet national safety standards
- **Traceable reprocessing cycles.** Document every cycle with: Date, patient name, procedure performed, serial number of RMD used, person reprocessing (connecting and removal), reprocessing method, cycle parameters (brand/batch number, expiry date, evidence of cycle result)

## TRANSDUCER REPROCESSING GUIDELINES

